



Safety Instructions

Read this safety instructions carefully before using the FLipMouse. Improper use might cause harm and injuries for which the AsTeRICS Foundation cannot be held responsible!

- **Do not** use this device if involuntary movements are possible (e.g.: spasms, epilepsy, Tourette syndrome, ...)
- **Do not share this device** with other users. If an infectious disease is indicated, replace the mouthpiece regularly (in accordance with a doctor).
- Clean the mouthpiece every 7 days ("Cleaning instructions" document)
- **Do not** use the device as mouth mouse if high salivation is indicated. Use the device as finger controlled mouse instead.
- **Do not** use the device in combination with safety critical devices or other medical devices.
- Use the device only with non-critical consumer electronics.
- **Do not** use the device in an explosive environment.